



# Laminar Air Flow Bench




**Model No. PLAFB-0124**

**Laminar Air Flow Bench** is a consistent tool provided with an HEPA - filtered work area which is essential to maintain dust particles and bacteria free clean air environment. The instrument is widely used in testing laboratories to maintain clean work area that provide protection and safety measures from outside harmful contaminants.


- Equipped with HEPA filter in work area
- Standardized and controlled lightening conditions
- Easy to use Menu driven system
- Wide Viewing Area
- High Quality Light sources used
- Plexi-glass frames for sample viewing




Top View



Front View



Side View



### Key Specifications:

Filter	HEPA filter
Lights	FL Fluorescent Light UV Blacklight Ultra Violet Black Light
Blower Motor	1/4 HP double shaft
Blower Timer	Optional (min / hrs)
Working area	4 x 2 feet (other size also available on request)

### Highlights:

Viewing Area Material	Transparent Plexiglass
No. of UV lights	1
No. of FL lights	1
Main Body Material	Steel
Dimensions	1200 x1005x 1515mm
Air Flow Gauge	Analogue Type

### Optional :

- Other sizes also available on request.

Thank you customers for choosing us as your partners in growth !



\*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

### OUR OTHER PRODUCTS :-



### Head Office : Presto Stantest Private Limited

I-42, DLF Industrial Area Phase-I, Delhi Mathura Road, Faridabad 121003, Haryana, India  
P : +91 129 4272727 E : info@prestogroup.com

### Branches :

#### Mumbai

503-504B, Fair link Centre, Opp. Fame Adlab  
Behind Mongini's Off Link Road  
Andheri (West), Mumbai 400058  
P : +91 022 42957052, 42957053  
E : west@prestogroup.com

#### Kolkata

93, 3<sup>rd</sup> Floor, Karnani Estate  
209, A.J.C. Bose Road  
Kolkata-700017. West Bengal  
P : +91 033 40612012  
E : east@prestogroup.com

#### Chennai

A-3, First Floor, Sri Sathya Sai Building  
755, Poonamallee High Road, Kilpauk  
Chennai- 600 010. Tamil Nadu  
P : +91 044 43865260, 43838372  
E : south@prestogroup.com

