






Laminar Air Flow Bench



Model No. PLAFB-0124

Laminar Air Flow Bench is a consistent tool provided with an HEPA - filtered work area which is essential to maintain dust particles and bacteria free clean air environment. The instrument is widely used in testing laboratories to maintain clean work area that provide protection and safety measures from outside harmful contaminants.

- Equipped with **HEPA filter** in work area
- **Standardized and controlled lightening conditions**
- Easy to use **Menu driven system**
- **Wide Viewing Area**
- High Quality Light sources used
- **Plexi-glass frames for sample viewing**





Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Sonipat • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad

www.prestogroup.com



Key Specifications:

| | |
|--------------|---|
| Filter | HEPA filter |
| Lights | FL Fluorescent Light UV Black light Ultra Violet Black Light |
| Blower Motor | 1/5 HP |
| Blower Timer | Hrs : min |

Highlights:

| | |
|-----------------------|--------------------------|
| Viewing Area Material | Transparent Plexi -glass |
| No. of UV lights | 1 |
| No. of FL lights | 1 |
| Main Body Material | Steel |
| Dimensions | 1200 x 1005 x 1515mm |

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-

| | | | | | |
|---|---|---|---|---|---|
|  |  |  |  |  |  |
| Top Load Tester | Wall Thickness Gauge | Gloss Meter Triangle | Eccentricity Tester | Moisture Analyzer MB 45 | Preform Neck Cutter |

