




Vibration Table



Model No. PVT-100

Vibration Table is a consistent tool for checking road worthiness of packages. It simulates the vibrations and jolts subjected to various packages and products during handling and transportation. PRESTO Vibration Table is manufactured under various Test Standards ASTM D999, TAPPI T-17, IS 7028(part II).

Available in Electro-mechanical, Pneumatic and Electromagnetic models also

- Digital LED Display
- Highly accurate test results under various frequency settings
- Digital Timer for setting test run time
- Sample Slippage protection by rail on all sides
- Strong base plate with rugged structure
- Adjustable Frequency
- Feather Touch Controls
- Short key feature available

The damage caused by repetitive movement of the objects is replicated by the Vibration Table. User can set frequency and time and examine the damage caused by visual inspection on completion of test.

